



Solutions for blistered medications

Aids for the storage and distribution of blistered medications in hospitals and care homes



Your expert for medication solutions

Wiegand has been developing convincing solutions for supplies and medication logistics for over 45 years. We offer tailor-made solutions - from central storage to the furnishing and equipping of ward rooms and pharmacies to the individual distribution of medicines at the patient's bedside. Wiegand does not offer standard solutions, but the right solution for every requirement.

Optimal storage and distribution of blistered medicines

The externally supplied and blistered medication significantly reduces the workload for the nursing staff. The challenge comes in the storage and delivery of the medication in blister rolls to the residents. Wiegand offers you a comprehensive range of aids to simplify these processes.

Wiegand's blister trolleys are tailored to the medication process with medication blisters and offer space for blister rolls as well as for additional medications such as creams, drops and supplementary materials.




Swiss quality products

As a family business, Wiegand will continue to rely on Switzerland as a production location in the future. The FlexLine trolleys are manufactured in the company's own production facility. Each FlexLine trolley is individually manufactured for you from high-quality and durable materials.

Easy storage, arrangement and distribution of blistered medicines

Wiegand offers a diverse range of aids for the storage and distribution of blistered medicines. Wiegand products simplify workflows and help to increase patient safety. The Wiegand system supports the distribution of blistered and other medicines in a simple, practical, and safe way.



therapies are changed or are not yet available in packet strip



Dispenser Marker

The Dispenser Marker is a permanent ink pen with special eraser. The fine tip is suitable for writing on medicine dispensers, medicine cups, medicine containers and labels.

WIEGAND® Starter-sets blistering

The starter sets support the distribution process of blistered pills in blister packes. The set includes one label holder per resident, a practical holder for the blister pack and a MediCup for dispensing the de-blistered pills. The necessary patient data and the patient picture can be printed on the enclosed labels and the labels can be stuck on the label holders provided. This increases medication safety, as an additional visual check can take place through the patient picture. The enclosed MediDispenser can be used for the medication of new patients.

Description Starter-set blistering S

The Starter Set Blistering S is ideal for the distribution of medication in blister pack for up to 5 residents.

The distribution tray has 10 grooves and 30 positions for MediCups.

Dimensions distribution tray: 260x300mm

Set consisting of:

1 distribution tray (S)

50 MediCups

100 lids for MediCups

5 holders for blister packs

5 Label holders with sliding lids 5 Sheets of label for patient picture and data

1 MediDispenser 341/349 with flexible subdivisions

Flyer: Informations for rest and nursing homes

Article No.



600.301.211.de 1 Starter-set blistering S

Description Starter-set blistering M

The Starter Set Blistering M is ideal for the distribution of medication in blister pack for up to 10 residents.

The distribution tray has 10 grooves and

50 positions for MediCups.

Dimensions distribution tray: 400x300mm

Set consisting of:

1 distribution tray (M)

50 MediCups, red

50 MediCups, yellow

100 lids for MediCups 10 holders for blister packs

10 Label holders with sliding lids

10 Sheets of label for patient picture and data

1 MediDispenser 341/349 with flexible subdivisions

1 Dispenser Marker, waterproof foil pen

1 MediSplitter WIEGAND®, tablet divider

Flyer: Informations for rest and nursing homes

Article No. 600.301.221.de 1 Starter-set blistering M







Description

Starter-set blistering L

The Starter Set Blistering L is ideal for the distribution of medication in blister pack for up to 16 residents. The distribution tray has 16 grooves and

80 positions for MediCups.

Dimensions distribution tray: 450x400mm

Set consisting of:

1 distribution tray (L)

50 MediCups, red

50 MediCups, yellow

100 lids for MediCups

16 holders for blister packs

16 Label holders with sliding lids

10 Sheets of label for patient picture and data

1 MediDispenser 341/349 with flexible subdivisions

1 Dispenser Marker, waterproof foil pen

1 MediSplitter WIEGAND®, tablet divider

Flyer: Informations for rest and nursing homes

Article No.

600.301.241.de

1 Starter-set blistering L







Description

Starter-set blistering XL

The Starter Set Blistering XL is ideal for the distribution of medication in blister pack for up to 22 residents.

The distribution tray has 22 grooves and

110 positions for MediCups.

Dimensions distribution tray: 400x600mm

Set consisting of:

1 distribution tray (XL)

50 MediCups, red

50 MediCups, yellow

100 lids for MediCups

22 holders for blister packs

22 Label holders with sliding lids 10 Sheets of label for patient picture and data

1 MediDispenser 341/349 with flexible subdivisions

1 PharmaBox small, small container with sliding lid

1 Dispenser Marker, waterproof foil pen

1 MediSplitter WIEGAND®, tablet divider

Flyer: Informations for rest and nursing homes



600.301.251.de


1 Starter-set blistering XL











Storage and distribution of medication in blister packs

Blistering is a variant of external preparation. This involves the mechanical preparation of prescribed tablets by a blister machine or device. The blister packs can be easily integrated into existing Wiegand distribution systems.

The hygienic and durable Wiegand BlisterBoxes store a patient's blister pack for one or two weeks and can be easily unrolled and detached.

Description BlisterBox duro

- Large dispenser for the storage and distribution of medicines blistered in a packet strip
- The BlisterBox duro is dishwasher safe, but can also be easily cleaned by hand
- Dimensions: 87x127x120 mm

Article No. 150.322.001 1 BlisterBox duro

150.322.001-54 Pack of 54 pcs. BlisterBox duro



Description BlisterBox Compact

- Large dispenser for the storage and distribution of medicines blistered in a packet strip
- The BlisterBoxes are delivered unfolded
- · Standard soap detergents can be used for cleaning
- Dimensions: 78x128x107 mm

Article No. 312T 1 BlisterBox Compact

312T-100 Pack of 100 pcs. BlisterBox Compact



Description BlisterBox Large

- Large dispenser for the storage and distribution of medicines blistered in a packet strip
- The BlisterBoxes are delivered unfolded
- · Standard soap detergents can be used for cleaning
- Dimensions: 78x162x139 mm

Article No. 311T 1 BlisterBox Large

311T-100 Pack of 100 pcs. BlisterBox Large



Wiegand distribution trays

Whether you're working with MediDispenser, blister packs or MediCup at the moment – the Wiegand distribution trays lend themselves to a wide range of uses. The low weight of the distribution tray and the flexible organisation support the medication dispensing process. Thanks to the clear overview this provides, administration errors are easily avoided



Description Distribution tray (S)

Stackable distribution tray with 10 grooves and 30 positions for MediCups

Suitable for Wiegand trolleys with trays

• Dimensions: 260x300 mm

Article No. 140.211 1 distribution tray (S)

140.211-45 Pack of 45 pcs. distribution trays (S)



Description Distribution tray (M)

- Stackable distribution tray with 10 grooves and 50 positions for MediCups
- Suitable for Wiegand trolleys with trays
- · Suitable for storage in ISO-A cabinets
- Dimensions: 400x300 mm

Article No. 140.221 1 distribution tray (M)

140.221-18 Pack of 18 pcs. distribution trays (M)



Easy storage and presentation of blistered medicines

Wiegand offers a diverse range of aids for the storage and distribution of blistered medicines. Wiegand products simplify workflows and help to increase patient safety. The Wiegand system enables the distribution of blistered and other medicines in a simple, practical and safe way.

Wiegand's hygienic and durable BlisterBoxes are high quality and robust and meet the high demands of nursing staff.

Description Holder for blister packs or labels

- Suitable for attaching tube blisters
- Pluggable on WIEGAND distribution trays,
- FlexModul trays/baskets
- Dimensions: 75x40 mm

Article No. 500.240XO.075 1 holder for blister bags or labels



Description Label holder with sliding lid

 Label holder as a carrier for patient data during medication dispensing with a Wiegand distribution tray


- Fits labels 510.13A4 and 510.13
- Dimensions: 77x55 mm

Art. No. 15/15SD-50 50 pc. label holder 15/15SD-500 500 pcs. label holder

Accessories 15SD.50 Spare sliding lids, 50 pcs.

510.13A4 Label for patient image and data

510.13 Label for medication



Description PharmaBox

- Ideal für the storage of blistered medicines and various auxiliary material
- With a universal holder (art. 16), it can easily be placed on a distribution tray
- Incl. sliding lid, transparent
- Dimensions: 134x72x30 mm

Art. No. 111/119.0-10 10 pcs. PharmaBox 111/119.0-400 400 pcs. PharmaBox


Accessories 119.40 Additional sliding lids, 40 pcs. 119.400 Additional sliding lids, 400 pcs.



Description Universal holder for distribution trays

- Universal holder for small containers and PharmaBoxes for mounting on Wiegand distribution trays
- Self-adhesive
- Dimensions: 30x65 mm

Art. No. 16-1 1 pc. universal holder 16-100 100 pcs. universal holder



Medication aids

Wiegand's wide and varied range offers everything needed for medication preparation, from pill processing aids to MediCups and storage boxes. Choose your suitable aids from our well-rounded range.

Description WIEGAND® PillCrusher Professional

- With this robust pill crusher, tablets can be processed hygienically and without the use of force to form a ready-to-use powder.
- Durable, easy to clean and maintenance-free
- Dimensions: 285x85x65 mm

Article No. 202.95.01 1 PillCrusher Professional

202.95.01-6 Pack of 6 pcs. PillCrusher Professional



Description

PillCrusher Pouches for single use

- The strong ribbed reinforced PillCrusher Pouches for single use prevent tearing, despite heavy pressures while crushing.
- With rounded corners so that the powder can be poured out easily
- Dimensions: 50x120 mm

Article No. 202.95.11-50 Pack of 50 pcs PillCrusher Pouches

202.95.11-1000 Pack of 1000 pcs. PillCrusher Pouches 202.95.11-8000 Pack of 8000 pcs. PillCrusher Pouches



Description

PillCrusher

- PillCrusher for the crushing and pulverisation of tablets
- Tablets can be processed directly in the medicine cup
- · Can be cleaned in the dishwasher
- Dimensions: 70x105 mm (Ø x H)

Article No. 202.097.001 1 Pillcrusher

202.097.001-10 Pack of 10 pcs. PillCrusher



Description

MediSplitter

- Handy tablet divider for halving oder quartering tablets
- The tablet can be easily positioned for division with the stainless steel blade
- Once split, the tablet falls into the storage container
- Dimensions: 85x35x25 mm

Article No. 96T 1 MediSplitter

96T-24 Pack of 24 pcs. MediSplitter



Description

Disposable MediCup ECO

- 40% less resource consumption
- With graduation up to 30 ml
- · Suitable for dispensing both oral medication and tablets
- For single use, not washable
- · Available in bars of 90 pcs. or in boxes of 9'000 pcs.
- Dimensions: 38x42 mm (Ø x H)

Art. No. 200.81.1 transparent 200.81.4 red 200.81.2 yellow 200.81.5 green

> 200.81.3 blue

Lid for MediCups ECO Accessories 200.82

Lid for MediCup ECO Description

• Lid fits on MediCup ECO

· Colour: white

• Dimensions: Ø 38 mm

Art. No. 200.82-100 100 pcs. lids 200.82-9000 9'000 pcs. lids



Disposable MediCup Light Description

- With graduation up to 30 ml
- · Suitable for dispensing both oral medication and tablets
- · For single use, not washable
- Available in bars of 80 pcs. or in boxes of 4'800 pcs.
- Dimensions: 38x42 mm (Ø x H)

Art. No. 200.77.1 transparent 200.77.4 red 200.77.2 vellow

> 200.77.3 blue

200.77.5 green

Lid for MediCups Light Accessories 200.78



Lid for disposable MediCup Light Description

· Lid fits on disposable MediCup Light

· Colour: white

• Dimensions: Ø 38 mm

Art. No. 200.78-100

100 pcs. lids 200.78-4000 4'000 pcs. lids



Description

MediCup reusable

- With graduation up to 25 ml
- · Suitable for dispensing both oral medication and tablets
- · Available in bars of 50 pcs. or in boxes of 4'000 pcs.
- Dimensions: 35x40 mm (Ø x H)

Art. No. 71.4 red transparent 71.2 yellow 71.5 green

> 71.3 blue

Accessories 72 Lid for MediCups reusable



Description

Lid for MediCup reusable

- · Writable lid fits on MediCup reusable
- Colour: transparent
- Dimensions: Ø 36 mm

Art. No. 72-100 100 pc. lids 72-2000 2'000 pcs. lids



Description

Dispenser for MediCups

- Dispenser for MediCups
- For mounting on wall or FlexLine trolleys
- · Without rod for attachment to trolley
- Material: Inox

Art. No. 64.746.081 1 pc. dispenser for MediCups ECO 64.746.077

1 pc. dispenser for MediCups Light 64.746.071 1 pc. dispenser for MediCups reusable



Description

Dispenser for lids for MediCups

- · Practical and hygienically storage container
- With Inox wall bracket
- · Separate opening for the direct removal of individual lids
- Dimensions: 195x185x234 mm

Art. No. 64.747 1 pc. dispenser



Description

Dish washer basket for MediCups reusable


- · Inox dish washer basket
- With space for 22 MediCups reusable
- or 18 inlay-dishes
- · Also available with two inserts
- Dimensions: 415x90x70 mm

Art. No. 910.71.1 910.66.1

1 pc. dish washer basket with insert for MediCups 1 pc. dish washer basket with insert for inlay-dishes

1 pc. dish washer basket with two inserts 910.71.66.1 910.00.1 1 pc. dish washer basket without inserts

910.07.1 1 pc. insert for MediCups 910.06.1 1 pc. insert for inlay-dishes



10

WIEGAND® blister trolleys

The Wiegand blister trolleys were configured together with specialists to suit the processes. The medication process with blistered medications requires certain working aids that make the process more efficient. The trolleys were equipped according to the needs of the professionals and can ideally support the work processes.

Wiegand's blister trolleys fulfil important criteria:

- Space-saving, hygienic and safe storage of medicines blistered in blister rolls
- Storage possibilities for additional medications and supplementary material
- Efficient, safe distribution and dispensing
- Time saving for the nursing staff, due to reduction of walking distances





Trolley for storage and distribution of card blisters

The Wiegand trolleys are also suitable for supporting medication processes that are implemented with blister cards. Thanks to the FlexLine system, it is possible to respond to individual needs and adapt the trolley.


Contact us for an individual consultation to configure the ideal trolley for your work process.

WIEGAND® FlexLine system


To tailor the blister trolley exactly to your needs, Wiegand offers you a high level of flexibility with the FlexLine system. Blister trolleys, like all other trolleys from Wiegand, can be individually configured thanks to our wide range of accessories and optimally adapted to your requirement.

Different variations available

- Wiegand blister trolleys are available in two widths: 45cm or 60cm
- Side panels or extraction frames
- Laterally extendable work surface
- · Additional shelves as storage space
- Various locking systems available
- Can be equipped with various accessories



Extraction frames for the fixed use of modules in the trolley



Side panels allow more flexible use of the trolley

Trolley colours

Whether you want to work with uniform colour codes or follow a specific furnishing concept - Wiegand offers you a diverse selection of trolley colours. The fronts of the drawers and doors are available in the following 10 standard colours at no extra charge:



The most popular blister trolleys

Blister trolley 20 residents

- 45 cm wide blister trolley with side panels
- Suitable for FlexModul trays 450x400x100 mm
- Space for up to six FlexModul trays
- Door made of hygienic and robust high pressure laminate
- Laterally extendable work surface
- Can be supplemented with the following accessories: Holder for MediCups dispenser, risk waste container, waste bag container, disinfectant dispenser and glove holder

Trolley-ID: B20.02002



Blister trolley 24 residents

- 45 cm wide blister trolley with extraction frames
- Suitable for FlexModul trays 450x400x100 mm
- Space for up to six FlexModul trays
- Door made of hygienic and robust high pressure laminate
- Laterally extendable work surface
- Can be supplemented with the following accessories: Holder for MediCups dispenser, risk waste container, waste bag container, disinfectant dispenser and glove holder

Trolley-ID: B20.02001



Blister trolley 30 residents

- 60 cm wide blister trolley with side panelsSuitable for FlexModul trays 600x400x100 mm
- Space for up to six FlexModul trays
- · Shelf for additional storage space or distribution tray shelf
- Door made of hygienic and robust high pressure laminate
- Laterally extendable work surface
- Can be supplemented with the following accessories: Holder for MediCups dispenser, risk waste container, waste bag container, disinfectant dispenser and glove holder

Trolley-ID: B20.02005



Blister trolley 30 residents

- 60 cm wide blister trolley with extraction framesSuitable for FlexModul trays 600x400x100 mm
- Shelf for additional storage space or distribution tray shelf
 Space up to five FlexModul trays
- Door made of hygienic and robust high pressure laminate
- Laterally extendable work surface
- Can be supplemented with the following accessories: Holder for MediCups dispenser, risk waste container, waste bag container, disinfectant dispenser and glove holder



Blister trolley 36 residents

- 60 cm wide blister trolley with extraction frames
- Suitable for FlexModul trays 600x400x100 mm
- Space for up to six FlexModul trays
- Door made of hygienic and robust high pressure laminate
- · Laterally extendable work surface
- Can be supplemented with the following accessories: Holder for MediCups dispenser, risk waste container, waste bag container, disinfectant dispenser and glove holder



FlexLine accessories

Wiegand offers you a selection of accessories which have been defined together with nursing professionals and ideally support the medication process with blistered medicines.

Description Lateral pull-out frame, with inserted tray

Additional work surface

Load up to 15 kg

Article No. 3440.2F Pull-out frame incl. tray



Description Holder with dispenser for disinfectant

· Also available without drainer

Article No. OP1416000 Disinfectant dispenser

64.263800 Drainer



Description Holder for hand gloves box Alu, adjustable

Practical glove holder for all popular glove packs

Article No. 64.711R-Alu Holder for hand gloves box right

64.711L-Alu Holder for hand gloves box left



Description Holder with risk waste container

• Risk waste containers are available in two sizes

Article No. 64.679-0.5 Risk waste container 0.5 Liter

64.679-2.0 Risk waste container 2 Liter

64.743 Holder for risk waste container 0.5 Litres 64.742 Holder for risk waste container 2 Litres



Description Waste bag holder

• with mounting / spring catch or fixed mounted available

Article No. 64.609F Waste bag holder large, swivelling (32x19x32cm)

64.610 Waste bag holder large, fix mounted (32x19x32cm)



Waste bag holder, Inox Description

· fix mounted

Article No. 64.752 Waste bag holder Inox (24x17x33cm)



FlexModul tray, light grey, Description

450x400x100 mm

· Openwork tray, divisible

Article No. 4540.210 Tray



Description

Subdivision set Acryl for FlexModul

tray 4540.210

· Ideal for storing directly in the trolley BlisterBoxes,

additional and reserve medications

4540.210-D2 Subdivision set Article No.



Description

FlexModul tray, light grey, 600x400x100 mm

· Openwork tray, divisible

Article No.

6040.210 Tray



Description

Subdivision set Acryl for FlexModul

tray 6040.210

 Ideal for storing directly in the trolley BlisterBoxes, additional and reserve

Article No. 6040.210-D2 Subdivision set



Description

Holder for blister bags or labels

Suitable for attaching tube blisters

· Pluggable on WIEGAND distribution trays, FlexModul

trays/baskets

• Dimensions: 75x40 mm

Article No.

500.240XO.075 1 holder for blister bags or labels



Locking systems for your Wiegand trolley

To safeguard your materials we off er various locking systems. Locking systems are usually chosen on the basis of a company-wide locking plan. We will be happy to help you choose a suitable system.

Description Preparation for digital half cylinder Euro

profile, 30 mm

· Central locking for trolley

· Preparation for on-site cylinder mounting

in drawer side panel or door

Article No. 45-1-L-VDigiZ

KitLock Description

· Electronic combination lock

KL1000* Article No.

Description **KABA 8 lock cylinder**

· Lock cylinder with bolt

· Bolt for front panel, door or side panel

Article No.

1031.20 (The precondition is a locking plan) **Alternative**









Cabinet for the storage of medicines and consumables

The storage of medicines and consumables in cabinets can - depending on the structural conditions - be used efficiently and offers various advantages for the organisation of the medication process. Our expert staff will accompany you on the way from planning to installation and will be at your side as a partner.





Efficient organisation ensures order in the medicine cabinets

Arrange your medicine cabinet in a space-saving and clear way. You can store the blister packet stips in the practical and hygienic BlisterBoxes and place the additional and reserve medication in FlexModul tays or TRESTON containers. Choose the ideal solution for your medication process.



TRESTON container

Description TRESTON container

• 1 front label incl. label protection included in the scope of delivery

Art. No. 3010-1 1 pc. TRESTON container, 92x300x82 mm

Accessories 10-D Divider

06-E-20 Label for divider, 20 pcs. 10-E-20 Front label, 20 pcs.

10-S-15 Protection for front label, 15 pcs.



Description TRESTON container

 1 front label incl. label protection included in the scope of delivery

Art. No. 3015-1 1 pc. TRESTON container, 132x300x100

Accessories 15-D Divider

06-E-20 Label for divider, 20 pcs. 15-E-20 Front label, 20 pcs.

15-S-15 Protection for front label, 15 pcs.



Description TRESTON container

• 1 front label incl. label protection included in the

scope of delivery

Art. No. 3020-1 1 pc. TRESTON container, 186x300x82 mm

Accessories 20-D Divider

06-E-20 Label for divider, 20 pcs.

20-E-20 Front label, 20 pcs.

20-S-15 Protection for front label, 15 pcs.



• 1 front label incl. label protection included in the

scope of delivery

Art. No. 4010-1 1 pc. TRESTON container, 92x400x82 mm

Accessories 10-D Divider

06-E-20 Label for divider, 20 pcs. 10-E-20 Front label, 20 pcs.

10-S-15 Protection for front label, 15 pcs.



Description TRESTON container

• 1 front label incl. label protection included in the

scope of delivery

Art. No. 4015-1 1 pc. TRESTON container, 132x400x100 mm

Accessories 15-D Divider

06-E-20 Label for divider, 20 pcs. 15-E-20 Front label, 20 pcs.

15-S-15 Protection for front label, 15 pcs.



Description TRESTON container

• 1 front label incl. label protection included in the

scope of delivery

Art. No. 4020-1 1 pc. TRESTON container, 186x400x82 mm

Accessories 20-D Divider

06-E-20 Label for divider, 20 pcs.

20-E-20 Front label, 20 pcs.

20-S-15 Protection for front label, 15 pcs.







Weekly trays or weekly containers for a single resident or storage trays for a whole group of patients

Wiegand offers you an extensive range of weekly systems which can be used individually. Store Medi-Dispensers practically on trays or in containers, depending on the space available. Wiegand solutions adapt to your work processes and offer many tools to improve patient safety.

Wiegand AG

Schlosserstrasse 5 CH-8180 Bülach Tel.: +41 44 872 74 74 wiegand@wiegand.ch www.wiegand.ch